

# Toilet Seat Wipes

Product # 1500

Make every public rest room toilet clean & fresh

With these wipes, you can be assured of a clean, healthy and useable toilet every time. Each wipe is premoistened with a potent formula that cleans on contact and evaporates within seconds. They come in 160-wipe canisters that fit into a dispenser that can be mounted anywhere.

## Use on:

- Toilet seats
- DO NOT** use as toilet tissue or for personal cleansing.

## Features and Benefits

- Fast-evaporating formula
- Quick and easy-to-use

## Wire Wipes Canister Holder - Small (product #9112)

- Quick and easy-to-use
- Use the wipes with a multi-purpose wipes dispenser, ideal as a toilet seat wipe.
- Wire, wall mount holds 77 and 81 mm canisters.
- The dispenser rapidly loads the 160-wipe canister.

## Properties

Appearance.....	White cloth
Fragrance.....	Alcohol
Size of towel.....	6" X 6"
Flash point (liquid).....	100 °F
VOC content (liquid).....	<10.0% by wt.

## Directions

Complete directions on product label  
Open canister feed the first towel through the opening. Insert canister upside down into dispenser.

## Ingredients

Ingredients	CAS Number
Water.....	7732-18-5
Isopropanol.....	67-63-0
N-Alkyl (68% C12, 32% C14) Dimethyl Ethyl Benzyl Ammonium Chloride.....	68956-79-6
N-Alkyl (60% C14, 30% C16, 5% C12, 5% C18) Dimethyl Benzyl Ammonium Chlorides.....	68391-01-5
Ethanol .....	64-17-5

## Safety


Keep out of reach of children.

## DOT Shipping (ground transportation)

Proper Shipping Name:	Consumer Commodity
Class:	ORM-D
ID Number:	None
Package Group:	None

## Quantities

160 wipes per canister                      6 per case



**TOTAL**  
SOLUTIONS®

Milwaukee, WI 53224 • 1-800-743-6417 • Fax 414-354-9219 • www.athea.com